

Multifunctional β-Ti₃AuAg Coatings For Enhanced Artificial Joints And Beyond

LEVERHULME TRUST_____


Cecil Cherian Lukose, Martin Birkett Email- c.c.lukose@northumbria.ac.uk




Joint Replacement Surgeries

Radically transforms the quality of life for patients suffering from degenerative joints, caused by various reasons But

It Is On The Rise



It Is Expensive



It Is Not Quick

Average Patient Waiting Time.

Year 2019 – 145

Year 2013 – 79

It Is Life Changing



- ➤ It costs NHS a good fortune £1.5 billion annually.
- > Physicians, surgeons and nursing teams are pushed beyond limit. > Waiting causes vulnerability & confidence issues.
- > Projected to be £2.4 billion by 2040, without accounting inflation. > Patient waiting time increase by 96.6% since 2019 level.
- Some categorise it as 'State Worse Than Death'.

Joint Replacement Surgery Is Best Done Once—But Joint Implants Do Fail !!!


Why Do Implants Fail.?

Ag doping

- > Almost 5-15% of joint implants fail within 10 years of primary -7500— 20,000 revisions/year.
- > Cost without inflation £75M/year to £200M/year Whereas complications cases cost even more.
- > Revision due to Prosthetic Joint Infection or PJI's on an average cost double the primary surgery.

How To Improve Their Lifecycle?

- > We coated articulating surface with multifunctional Ti₃Au-Ag thin films.
- > Ti₃Au formed superhard-wear resistant- biocompatible protection for these surfaces and Silver (Ag) provided protection against post surgery infections.




Top reasons for early hip and knee revisions

Prosthetic Joint infection

Peri-Prosthetic Fracture

Instability



Implant coating



How Can Implants Be Coated And Tested


Year 2021 – 284

> Articulating surface of artificial joints can be coated using Magnetron Sputtering Magnetron Sputtering technique. Ti6Al4V substrate

Colorimetric reactior


➤ Ti₃Au-Ag thin films were deposited on glass an Ti64 coupons.

> We explored the effects of varying -



Antimicrobial properties Mechanical properties Mechanical characterization ₹ 72 and 168 hours immersed in DMEM

4. Reagent Exposur




 Ar^+ Ar^+


Sputtering bias

The Performance Of Ti₃Au-Ag Multifunctional Coatings

Increased surface Hardness > 14 GPa




60 times increase in wear resistance



Wear Resistance

Cell viability greater than 95%



Biocompatibility

og Reduction of bacterial (a.u) Copper control Time (min)

Bacterial reduction >1.5log


Antimicrobial Resistance

Foundational Study

Surface Hardness



Clinical Trials







Publication list!

